

Pro-Care[®] CraigBed (Modular Rehab Enclosure Bed, Wall Bed, Cubicle Bed)



CraigBed with 4' Riser Block



CraigBed - Closed and Latched



CraigBed - Floor Pad

Modular Rehab Treatment Bed

The Pro-Care[®] CraigBed is a functional, reliable solution for facilitating the safe treatment of TBI - Traumatic Brain Injury, Huntington Disease, Cerebral Palsy, and other conditions. <u>As an alternative to using strap restraints</u>, the occupant may be placed in this padded environment with reduced likelihood of injury. It provides a safer solution to fully netted bed enclosures that carry the inherent risk of the occupant becoming entangled in the netting. The unit is composed of a foam mattress and four padded walls that connect to form a safe padded cubicle enclosure. The reversible, heavy duty foam mattress is covered in a tough, nylon reinforced hospital grade vinyl fabric. It is fluid proof, stain & static resistant, antimicrobial, durable, and easy to clean. The mattress may be placed on an optional wood platform, foam riser block or floor pad, or it may be placed directly on the floor.* The walls connect with web straps and quick-release buckles. The interior and exterior are covered in hospital grade vinyl that wipes clean with standard cleaning solutions.

Quick & Easy Access

Access to the interior of the bed is quick and easy from either of two opposite corners. Simply squeeze the quick release buckles and the walls swing open. To close again, loosen the adjustable straps to allow some slack, then snap the buckles back together and cinch up till the walls are again snug against the mattress. All hardware is reversible so that the entry can be changed from left to right side to facilitate your room layout.

Easy Handling

The Pro-Care[®] CraigBed is designed for easy transport and storage. Each wall section rolls on heavy duty, non-marring, lockable swivel casters. The two wall sections may be latched together for ease of handling.

Options: Riser Blocks, Floor Pads

Vinyl covered floor are available in 1" to 4" thicknesses to provide a barrier between the mattress and the floor without raising the mattress height significantly. A special platform is available for use with Hoyer Lifts.

Note!

The use of a riser block or platform **is NOT recommended** for active occupants over 6 feet tall, due to the possible potential for them to accidentally tumble over the top of the restraining walls when standing - even when using the Tall Walled Bed.

Flammability

Mattress meets Federal Open-Flame Flammability Standard CPSC 16/CFR/Part1633, and Cigarette Standard Part 1632. All materials used in the manufacture of this mattress meet the California Open Flame Standard T.B. #117.

* Check with local health and safety codes to see if the mattress may be placed directly on the floor.

v06.2025



Pro-Care[®] CraigBed (Modular Rehab Enclosure Bed, Wall Bed, Cubicle Bed)

	Standard-XL	Standard	Twin-XL	Twin
Outside				
Length	88	82	88	82
Width	56	56	45	45
Standard Walls CR420 Height	451/2	451/2	451/2	451/2
Tall Walls CR480 Height	511/2	511/2	511/2	511/2
Inside				
Length	80	74	80	74
Width	48	48	37	37
Standard Walls CR420 Height	391/2	391/2	391/2	391/2
Tall Walls CR480 Height	451/2	451/2	451/2	451/2
Mattress				
Length	80	74	80	74
Width	48	48	37	37
Thickness	6	6	6	6
Decubitus Mattress Thickness	51/2	51/2	51/2	51/2
		caster. The insi	red from the floo de wall height is r nattress.	

Order Information		
Includes 4 padded walls and a mattress. (Two Double-Walled hinged sections)		
Standard-XL size bed complete with 4 wall enclosure and 6" reversible mattress		
Standard size bed complete with 4 wall enclosure and 6" reversible mattress		
Twin-XL size bed complete with 4 wall enclosure and 6" reversible mattress		
Twin size bed complete with 4 wall enclosure and 6" reversible mattress		
REQUIRED when a riser block or platform is used. Adds 6" to wall height.		
Standard-XL size bed with Tall Walls complete with 4 wall enclosure and 6" reversible mattress		
Standard size bed with Tall Walls complete with 4 wall enclosure and 6" reversible mattress		
Twin-XL size bed with Tall Walls complete with 4 wall enclosure and 6" reversible mattress		
Twin size bed with Tall Walls complete with 4 wall enclosure and 6" reversible mattress		
When ordering a CraigBed with the 5 ¹ / ₂ " Decubitus Reducing Mattress, specify Model CR424-xxxx, or CR484-xxxx.		
Zippered fluid resistant covers are available to fit the CraigBed Mattresses. Item - FRZCN		

v06.2025



Pro-Care[®] CraigBed (Modular Rehab Enclosure Bed, Wall Bed, Cubicle Bed)

Order Information			
4" Foam Riser Block**	Optional vinyl covered foam riser block - raises mattress 4" for better viewing.		
CR124-4880	Standard-XL size 4" riser block		
CR124-4874	Standard size 4" riser block		
CR124-3780	Twin-XL size 4" riser block		
CR124-3774	Twin size 4" riser block		
1'' Floor Pad	Optional 1" vinyl covered foam floor pad to prevent mattress from lying directly on the floor.		
CR110-4880	Standard-XL size 4" 1" floor pad		
CR110-4874	Standard size 1" floor pad		
CR110-3780	Twin-XL size 1" floor pad		
CR110-3774	Twin size 1" floor pad		
Wood Platform**	Commercial platform to support the mattress.		
CR135-4880	Platform for use with Standard-XL size CraigBed		
CR135-4874	Platform for use with Standard size CraigBed		
CR135-3780	Platform for use with Twin-XL size CraigBed		
CR135-3774	Platform for use with Twin size CraigBed		
Hoyer Lift Platform**	Commercial platform for use with Hoyer Lift		
CR145-4880	Hoyer Lift Platform for use with Standard-XL size CraigBed		
CR145-4874	Hoyer Lift Platform for use with Standard size CraigBed		
CR145-3780	Hoyer Lift Platform for use with Twin-XL size CraigBed		
CR145-3774	Hoyer Lift Platform for use with Twin size CraigBed		
Zippered Covers	Fluid resistant nylon – Heavy Duty Zipper, Fire Retardant		
FRZCN-4880	Zippered Mattress Cover for use with Standard-XL size mattress		
FRZCN-4874	Zippered Mattress Cover for use with Standard size mattress		
FRZCN-3780	Zippered Mattress for use with Twin-XL size mattress		
FRZCN-3774	Zippered Mattress for use with Twin size mattress		
	**Requires Tall Walled CraigBed (models 480 or 484). Unsafe for use with standard height walls!		

v06.2025